
Summary

Production Name	Flotillin 1 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	FLOT1
Alternative Names	Flotillin-1
Gene ID	10211
SwissProt ID	O75955

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 47 kDa; Observed MW: 47 kDa

Background

Product Name: Flotillin 1 Rabbit Monoclonal Antibody
Catalog #: AMRe03002




May act as a scaffolding protein within caveolar membranes, functionally participating in formation of caveolae or caveolae-like vesicles.

Research Area

Cell Biology

Image Data



Western blot analysis of Flotillin 1 in HeLa, A549, MCF-7, U251 lysates using Flotillin 1 antibody.

Note

For research use only.