
Summary

Production Name	HDLBP Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	HDLBP
Alternative Names	HBP; VGL; PRO2900
Gene ID	3069
SwissProt ID	Q00341

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 141 kDa; Observed MW: 141 kDa

Background

Product Name: HDLBP Rabbit Monoclonal Antibody
Catalog #: AMRe03021




Appears to play a role in cell sterol metabolism.

Research Area

Cardiovascular

Image Data



Western blot analysis of HDLBP in HeLa, A549, HL-60, U251, U87-MG lysates using HDLBP antibody.

Note

For research use only.