

Product Name: LINE 1 ORF1p Rabbit Monoclonal Antibody
Catalog #: AMRe03068

Summary

Production Name	LINE 1 ORF1p Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	L1RE1
Alternative Names	L1ORF1p; L1RE1; LRE1
Gene ID	105370531.0
SwissProt ID	Q9UN81

Application

Dilution Ratio	WB: 1/500-1/1000 IP: 1/20
Molecular Weight	Calculated MW: 40 kDa; Observed MW: 40 kDa

Product Name: LINE 1 ORF1p Rabbit Monoclonal Antibody
Catalog #: AMRe03068


Background

Nucleic acid-binding protein which is essential for retrotransposition of LINE-1 elements in the genome.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of LINE1 ORF1p in HeLa, A549 lysates using LINE 1 ORF1p antibody.

Note

For research use only.