

Product Name: MPP1 Rabbit Monoclonal Antibody
Catalog #: AMRe03092



Summary

Production Name	MPP1 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	MPP1
Alternative Names	MRG1; PEMP; AAG12; EMP55; DXS552E
Gene ID	4354
SwissProt ID	Q00013

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 52 kDa; Observed MW: 55 kDa

Background

Product Name: MPP1 Rabbit Monoclonal Antibody
Catalog #: AMRe03092




Essential regulator of neutrophil polarity.

Research Area

-

Image Data



Western blot analysis of MPP1 in A549, HL-60 lysates using MPP1 antibody.

Note

For research use only.