
Summary

Production Name	Nectin 2 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,ICC/IF,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	NECTIN2
Alternative Names	HVEB; PRR2; CD112; PVRR2
Gene ID	5819
SwissProt ID	Q92692

Application

Dilution Ratio	WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20
Molecular Weight	Calculated MW: 58 kDa; Observed MW: 70-80 kDa

Background

Product Name: Nectin 2 Rabbit Monoclonal Antibody
Catalog #: AMRe03105




Probable cell adhesion protein.


Research Area

Immunology

Image Data



Western blot analysis of Nectin 2 in HeLa, A549, HL-60, U251, U87-MG lysates using Nectin 2 antibody.



Immunocytochemistry analysis of Nectin 2 (green) in heLa using Nectin 2 antibody, and DAPI (blue)

Note

For research use only.