
Summary

Production Name	RGS6 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	RGS6
Alternative Names	regulator of G protein signaling 6; GAP; S914; HA117
Gene ID	9628
SwissProt ID	P49758

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 54 kDa; Observed MW: 54 kDa

Background

Product Name: RGS6 Rabbit Monoclonal Antibody
Catalog #: AMRe03185




Regulates G protein-coupled receptor signaling cascades.

Research Area

-

Image Data



Western blot analysis of RGS6 in A549, HL-60, U2OS lysates using RGS6 antibody.

Note

For research use only.