

Product Name: DAP kinase 3 Rabbit Monoclonal Antibody
Catalog #: AMRe03261

Summary

Production Name	DAP kinase 3 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	DAPK3
Alternative Names	DAPK3; ZIPK; Death-associated protein kinase 3; DAP kinase 3; DAP-like kinase; Dlk; MYPT1 kinase; Zipper-interacting protein kinase; ZIP-kinase
Gene ID	1613
SwissProt ID	O43293

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 53 kDa; Observed MW: 53 kDa

Product Name: DAP kinase 3 Rabbit Monoclonal Antibody
Catalog #: AMRe03261


Background

Death-associated protein kinase 3 (DAPK3) induces morphological changes in apoptosis when overexpressed in mammalian cells. These results suggest that DAPK3 may play a role in the induction of apoptosis.

Research Area

Cell Biology

Image Data



Western blot analysis of ZIP Kinase in HeLa, A549, MCF-7, U251 lysates using DAP kinase 3 antibody.

Note

For research use only.