

**Product Name: MonoMethyl-Histone H3 (Lys18) Rabbit
Monoclonal Antibody
Catalog #: AMRe03281**

Summary

Production Name	MonoMethyl-Histone H3 (Lys18) Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Methylated
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	H3C1
Alternative Names	H3K18me; H3/j; H3C1; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FJ; H3C10; H3C11; HIST1H3J
Gene ID	8350
SwissProt ID	P68431

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 15 kDa; Observed MW: 15 kDa

**Product Name: MonoMethyl-Histone H3 (Lys18) Rabbit
Monoclonal Antibody
Catalog #: AMRe03281**


Background

H3 Core component of nucleosome. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability.


Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of MonoMethyl-Histone H3 (Lys18) in 3T3, rat brain lysates using MonoMethyl-Histone H3 (Lys18) antibody.



Western blot analysis of MonoMethyl-Histone H3 (Lys18) in HeLa, A549, HL-60, U2OS, U87-MG lysates using MonoMethyl-Histone H3 (Lys18) antibody.

Note

For research use only.