

Summary

Production Name	PRPF31 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,ICC/IF
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	PRPF31
Alternative Names	Protein 61K; hPrp31; U4/U6 snRNP 61 kDa protein
Gene ID	26121
SwissProt ID	Q8WWY3

Application

Dilution Ratio	WB: 1/500-1/1000 IF: 1/50-1/200
Molecular Weight	Calculated MW: 55 kDa; Observed MW: 55 kDa

Background

Product Name: PRPF31 Rabbit Monoclonal Antibody
Catalog #: AMRe03282




Involved in pre-mRNA splicing as component of the spliceosome.


Research Area

Epigenetics and Nuclear Signaling


Image Data



Western blot analysis of PRPF31 in PC-12 lysates using PRPF31 antibody.



Western blot analysis of PRPF31 in C2C12, PC-3 lysates using PRPF31 antibody.



Immunocytochemistry analysis of PRPF31 (green) in HeLa using PRPF31 antibody, and DAPI (blue)

Product Name: PRPF31 Rabbit Monoclonal Antibody
Catalog #: AMRe03282



Note

For research use only.