

**Product Name: Phospholipase A2 Activator Protein
Rabbit Monoclonal Antibody
Catalog #: AMRe03291**

Summary

Production Name	Phospholipase A2 Activator Protein Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,ICC/IF
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	PLAA
Alternative Names	Glycerophosphatase; PLA2P; Plaa; PLAP
Gene ID	9373
SwissProt ID	Q9Y263

Application

Dilution Ratio	WB: 1/500-1/1000 IF: 1/50-1/200
Molecular Weight	Calculated MW: 87 kDa; Observed MW: 87 kDa

**Product Name: Phospholipase A2 Activator Protein
Rabbit Monoclonal Antibody
Catalog #: AMRe03291**


Background

Involved in the maintenance of ubiquitin levels.


Research Area

Immunology

Image Data



Immunocytochemistry analysis of Phospho-lipase A2 Activator Protein (green) in HeLa using Phospho-lipase A2 Activator Protein antibody, and DAPI (blue).



Western blot analysis of Phospho-lipase A2 activator protein (PLAP) in 3T3, HeLa lysates using Phospho-lipase A2 activator protein (PLAP) antibody.

Note

For research use only.