

**Product Name: Adiponectin Rabbit Monoclonal Antibody**  
**Catalog #: AMRe03758**

---

## Summary

<b>Production Name</b>	Adiponectin Rabbit Monoclonal Antibody
<b>Description</b>	Recombinant Rabbit Monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC-F,IHC-P,ICC/IF
<b>Reactivity</b>	Mouse,Rat

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal Antibody
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
<b>Purification</b>	Affinity Purified

## Immunogen

<b>Gene Name</b>	Adipoq
<b>Alternative Names</b>	Ad; APN; Acdc; Adid; apM1; 30kDa; GBP28; adipo; Acrp30
<b>Gene ID</b>	11450.0
<b>SwissProt ID</b>	Q60994

## Application

<b>Dilution Ratio</b>	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200
<b>Molecular Weight</b>	Calculated MW: 27 kDa; Observed MW: 30 kDa

**Product Name: Adiponectin Rabbit Monoclonal Antibody**  
**Catalog #: AMRe03758**


## Background

Important adipokine involved in the control of fat metabolism and insulin sensitivity, with direct anti-diabetic, anti-atherogenic and anti-inflammatory activities. Stimulates AMPK phosphorylation and activation in the liver and the skeletal muscle, enhancing glucose utilization and fatty-acid combustion. Antagonizes TNF-alpha by negatively regulating its expression in various tissues such as liver and macrophages, and also by counteracting its effects. Inhibits endothelial NF-kappa-B signaling through a cAMP-dependent pathway. May play a role in cell growth, angiogenesis and tissue remodeling by binding and sequestering various growth factors with distinct binding affinities, depending on the type of complex, LMW, MMW or HMW.


## Research Area

Cardiovascular

## Image Data



Western blot analysis of Adiponectin in rat fat lysates using Adiponectin antibody.



Western blot analysis of Adiponectin in 3T3-L1 lysates using Adiponectin antibody.

## Note

For research use only.

**Product Name: Adiponectin Rabbit Monoclonal  
Antibody**  
**Catalog #: AMRe03758**

---

