

Product Name: YTHDF1 Rabbit Monoclonal Antibody
Catalog #: AMRe03818



Summary

Production Name	YTHDF1 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Hamster

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	Ythdf1
Alternative Names	C20orf21; DF1
Gene ID	228994.0
SwissProt ID	P59326

Application

Dilution Ratio	WB: 1/500-1/1000 IP: 1/20
Molecular Weight	Calculated MW: 61 kDa; Observed MW: 61 kDa

Background

Product Name: YTHDF1 Rabbit Monoclonal Antibody
Catalog #: AMRe03818




Specifically recognizes and binds N6-methyladenosine (m6A)-containing mRNAs, and regulates their stability.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of YTHDF1 in HeLa, CHO-K1, Ramos lysates using YTHDF1 antibody.

Note

For research use only.