

Summary

Production Name	SESN1 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,FC
Reactivity	Human, Mouse, Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.
Purification	Affinity Purified

Immunogen

Gene Name	SESN1
Alternative Names	PA26; sesn1; SEST1; sestrin 1
Gene ID	27244
SwissProt ID	Q9Y6P5

Application

Dilution Ratio	WB: 1/500-1/1000 FC: 1/50-1/100
Molecular Weight	Calculated MW:57 kDa;Observed MW: 70 kDa

Background

Product Name: SESN1 Rabbit Monoclonal Antibody
Catalog #: AMRe03914




Involved in the reduction of peroxiredoxins. May also be regulator of cellular growth.

Research Area

-

Image Data



Western blot analysis of SESN1 in HeLa,K562 lysates using SESN1 antibody.

Note

For research use only.