

Summary

Production Name	PPIC Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.
Purification	Affinity Purified

Immunogen

Gene Name	PPIC
Alternative Names	PPIC,CYPC,Rotamase C,Cyclophilin C,Peptidyl-prolyl cis-trans isomerase C
Gene ID	5480.0
SwissProt ID	P45877

Application

Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100
Molecular Weight	Calculated MW:23 kDa;Observed MW: 23 kDa

Background

Product Name: PPIC Rabbit Monoclonal Antibody
Catalog #: AMRe03949




PPIase that catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides and may therefore assist protein folding.

Research Area

Immunology

Image Data



Western blot analysis of PPIC in HCT116 lysates using PPIC antibody.

Note

For research use only.