

Summary

Production Name	Syntaxin 3 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P,ICC/IF
Reactivity	Human, Mouse, Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.
Purification	Affinity Purified

Immunogen

Gene Name	STX3
Alternative Names	Stx3; STX3A; Syntaxin 3A
Gene ID	6809
SwissProt ID	Q13277

Application

Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200
Molecular Weight	Calculated MW:33 kDa;Observed MW: 33 kDa

Background

Product Name: Syntaxin 3 Rabbit Monoclonal Antibody
Catalog #: AMRe03964




Potentially involved in docking of synaptic vesicles at presynaptic active zones.

Research Area

Neuroscience

Image Data



Western blot analysis of Syntaxin 3 in HepG2, Raw 264.7 lysates using Syntaxin 3 antibody.

Note

For research use only.