

Product Name: PI 3 Kinase p150 Rabbit Monoclonal Antibody
Catalog #: AMRe03969



Summary

Production Name	PI 3 Kinase p150 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,ICC/IF
Reactivity	Human, Mouse, Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.
Purification	Affinity Purified

Immunogen

Gene Name	PIK3R4
Alternative Names	P150; Phosphoinositide 3 kinase regulatory subunit 4; PI3 kinase p150 subunit; PIK3R4; VPS15
Gene ID	30849
SwissProt ID	Q99570

Application

Dilution Ratio	WB: 1/500-1/1000 IF: 1/50-1/200
Molecular Weight	Calculated MW:153 kDa;Observed MW: 153 kDa

Product Name: PI 3 Kinase p150 Rabbit Monoclonal Antibody
Catalog #: AMRe03969




Background

Regulatory subunit of the PI3K complex. May regulate membrane trafficking late in the endocytic pathway.

Research Area

Signal Transduction

Image Data



Western blot analysis of PI 3 Kinase p150 in Jurkat,Raji lysates using PI 3 Kinase p150 antibody.

Note

For research use only.