

Product Name: Sulfotransferase 2A1 Rabbit Monoclonal Antibody
Catalog #: AMRe03998

Summary

Production Name	Sulfotransferase 2A1 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P,ICC/IF,FC
Reactivity	Human, Mouse, Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.
Purification	Affinity Purified

Immunogen

Gene Name	SULT2A1
Alternative Names	HST; ST2; STD; hSTa; DHEAS; ST2A1; ST2A3; DHEA-ST
Gene ID	6822
SwissProt ID	Q06520

Application

Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 FC: 1/50-1/100
Molecular Weight	Calculated MW:34 kDa;Observed MW: 34 kDa

Product Name: Sulfotransferase 2A1 Rabbit Monoclonal Antibody
Catalog #: AMRe03998




Background

Catalyzes the sulfation of steroids and bile acids in the liver and adrenal glands.

Research Area

Signal Transduction

Image Data



Western blot analysis of Sulfotransferase 2A1 in Human fetal liver lysates using Sulfotransferase 2A1 antibody.

Note

For research use only.