

Summary

Production Name	CD53 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P,FC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.
Purification	Affinity Purified

Immunogen

Gene Name	CD53
Alternative Names	CD53; MOX44; TSPAN25; Leukocyte surface antigen CD53; Cell surface glycoprotein CD53; Tetraspanin-25; Tspan-25; CD53
Gene ID	963
SwissProt ID	P19397

Application

Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100 FC: 1/50-1/100
Molecular Weight	Calculated MW:24 kDa;Observed MW: 35-55kDa

Product Name: CD53 Rabbit Monoclonal Antibody
Catalog #: AMRe04003




Background

May be involved in growth regulation in hematopoietic cells.

Research Area

Immunology

Image Data



Western blot analysis of CD53 in Raji lysates using CD53 antibody.

Note

For research use only.