

Product Name: alpha Actinin/ACTN1 Rabbit Monoclonal Antibody
Catalog #: AMRe04010

Summary

Production Name	alpha Actinin/ACTN1 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human, Mouse, Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.
Purification	Affinity Purified

Immunogen

Gene Name	ACTN1
Alternative Names	actinin 1 smooth muscle; actinin alpha 1; ACTN 1; FLJ40884
Gene ID	87
SwissProt ID	P12814

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW:103 kDa;Observed MW: 103 kDa

Product Name: alpha Actinin/ACTN1 Rabbit Monoclonal Antibody
Catalog #: AMRe04010




Background

F-actin cross-linking protein which is thought to anchor actin to a variety of intracellular structures. This is a bundling protein.

Research Area

Signal Transduction

Image Data



Western blot analysis of alpha Actinin/ACTN1 in 3T3 lysates using alpha Actinin/ACTN1 antibody.

Note

For research use only.