

---

## Summary

<b>Production Name</b>	Caspase 14 Rabbit Monoclonal Antibody
<b>Description</b>	Recombinant Rabbit Monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC-F,IHC-P,ICC/IF,FC,IP
<b>Reactivity</b>	Human, Mouse, Rat

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal Antibody
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.
<b>Purification</b>	Affinity Purified

## Immunogen

<b>Gene Name</b>	CASP14
<b>Alternative Names</b>	CASP14; Caspase-14; Caspase 14; MICE; CASP-14
<b>Gene ID</b>	23581
<b>SwissProt ID</b>	P31944

## Application

<b>Dilution Ratio</b>	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20 FC: 1/50-1/100
<b>Molecular Weight</b>	Calculated MW:28 kDa;Observed MW: 28 kDa

**Product Name: Caspase 14 Rabbit Monoclonal Antibody**  
**Catalog #: AMRe04015**




## Background

Caspases are a family of cysteine proteases that play an essential role in carrying out apoptosis. Caspase-14, also named MICE, is a unique member of the caspase family with restricted expression; it is found in embryonic tissues and adult skin. Caspase-14 is weakly processed into p18 and p11 subunits by caspase-8. May also be responsible for proteolytic processing of filaggrin during terminal differentiation of keratinocytes.

## Research Area

Cell Biology

## Image Data



Western blot analysis of Caspase 14 in MCF7, HACAT lysates using Caspase 14 antibody.

## Note

For research use only.