
Summary

Production Name	GSTA4 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human, Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.
Purification	Affinity Purified

Immunogen

Gene Name	GSTA4
Alternative Names	GSTA4;
Gene ID	2941.0
SwissProt ID	O15217

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW:26 kDa;Observed MW: 26 kDa

Background

Product Name: GSTA4 Rabbit Monoclonal Antibody
Catalog #: AMRe04019




Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles. This isozyme has a high catalytic efficiency with 4-hydroxyalkenals such as 4-hydroxynonenal (4-HNE).

Research Area

Signal Transduction

Image Data



Western blot analysis of GSTA4 in 293T lysates using GSTA4 antibody.

Note

For research use only.