

Product Name: Dishevelled 3 Rabbit Monoclonal Antibody
Catalog #: AMRe04028

Summary

Production Name	Dishevelled 3 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P,ICC/IF,FC
Reactivity	Human, Mouse, Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.
Purification	Affinity Purified

Immunogen

Gene Name	DVL3
Alternative Names	DRS3,KIAA0208,Dishevelled-3,DSH homolog 3,Segment polarity protein dishevelled homolog DVL-3
Gene ID	1857.0
SwissProt ID	Q92997

Application

Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 FC: 1/50-1/100
Molecular Weight	Calculated MW:78 kDa;Observed MW: 95 kDa

Product Name: Dishevelled 3 Rabbit Monoclonal Antibody
Catalog #: AMRe04028


Background

Involved in the signal transduction pathway mediated by multiple Wnt genes.

Research Area

Neuroscience

Image Data



Western blot analysis of Dishevelled 3 in K562 lysates using Dishevelled 3 antibody.

Note

For research use only.