

Product Name: CRABP2 Rabbit Monoclonal Antibody
Catalog #: AMRe04040



Summary

Production Name	CRABP2 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,ICC/IF
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	CRABP2
Alternative Names	Cellular retinoic acid-binding protein II; CRABP-II
Gene ID	1382
SwissProt ID	P29373

Application

Dilution Ratio	WB: 1/500-1/1000 IF: 1/50-1/200
Molecular Weight	Calculated MW: 16 kDa; Observed MW: 16 kDa

Background

Product Name: CRABP2 Rabbit Monoclonal Antibody
Catalog #: AMRe04040




Transports retinoic acid to the nucleus.


Research Area

Signal Transduction

Image Data



Western blot analysis of CRABP2 in HeLa lysates using CRABP2 antibody.



Immunocytochemistry analysis of CRABP2 (green) in MCF-7 using CRABP2 antibody, and DAPI (blue)

Note

For research use only.