

**Product Name: Acetyl Coenzyme A Carboxylase (6N12)
Rabbit Monoclonal Antibody
Catalog #: AMRe04158**



Summary

Production Name	Acetyl Coenzyme A Carboxylase (6N12) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Acetyl Antibody
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% New type preservative N and 0.05% BSA
Purification	Affinity purification

Immunogen

Gene Name	ACACB
Alternative Names	
Gene ID	32.0
SwissProt ID	O00763.

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	277 kDa


Background

**Product Name: Acetyl Coenzyme A Carboxylase (6N12)
Rabbit Monoclonal Antibody
Catalog #: AMRe04158**



Research Area

Image Data



Western blot analysis of Acetyl Coenzyme A Carboxylase in HeLa, HCT116 lysates using Acetyl Coenzyme A Carboxylase antibody.

Note

For research use only.