

Product Name: Acetyl-Histone H3 (Lys27) (17F16)
Rabbit Monoclonal Antibody
Catalog #: AMRe04172

Summary

Production Name	Acetyl-Histone H3 (Lys27) (17F16) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Acetyl Antibody
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	H3C1
Alternative Names	
Gene ID	8350.0
SwissProt ID	P68431.A synthetic acetyl-peptide corresponding to residues surrounding Lys27 of human Histone H3

Application


Dilution Ratio	WB: 1:1000
Molecular Weight	15kDa

Product Name: Acetyl-Histone H3 (Lys27) (17F16)
Rabbit Monoclonal Antibody
Catalog #: AMRe04172

Background

Research Area

Image Data



Western blot detection of Histone H3 (acetyl-Lys27) in HeLa,A549,HL-60 cell lysates using Histone H3 (acetyl-Lys27) antibody(1:1000 diluted).

Note

For research use only.