

**Product Name: Phospho-CDK1/2 (Thr14) (5Z8) Rabbit
Monoclonal Antibody
Catalog #: AMRe05871**



Summary

Production Name	Phospho-CDK1/2 (Thr14) (5Z8) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Phospho Antibody
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	CDK1/CDK2/
Alternative Names	
Gene ID	983/1017
SwissProt ID	P06493/P24941.A synthetic phosphopeptide corresponding to residues surrounding Thr14 of human Cdk1/2

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	34kDa


**Product Name: Phospho-CDK1/2 (Thr14) (5Z8) Rabbit
Monoclonal Antibody
Catalog #: AMRe05871**



Background

Research Area

Image Data



Western blot detection of Cdk1/2 (phospho-Thr14) in HepG2,U2OS,U251 using Cdk1/2 (phospho-Thr14) antibody(1:1000 diluted)

Note

For research use only.