

**Product Name: Phospho-eIF4B (Ser406) (13H6) Rabbit
Monoclonal Antibody
Catalog #: AMRe05893**



Summary

Production Name	Phospho-eIF4B (Ser406) (13H6) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Phospho Antibody
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	EIF4B
Alternative Names	
Gene ID	1975.0
SwissProt ID	P23588.A synthetic phosphopeptide corresponding to residues surrounding Ser406 of human eIF4B

Application

Dilution Ratio	WB: 1:2000-1:10000
Molecular Weight	69kDa


**Product Name: Phospho-eIF4B (Ser406) (13H6) Rabbit
Monoclonal Antibody
Catalog #: AMRe05893**



Background

Research Area

Image Data



Western blot detection of eIF4B (phospho-Ser406) in A549,HL-60,U2OS,C6 using eIF4B (phospho-Ser406) antibody(1:1000 diluted)

Note

For research use only.