

**Product Name: Phospho-FAK (Tyr576) (17X9) Rabbit  
Monoclonal Antibody  
Catalog #: AMRe05905**



## Summary

<b>Production Name</b>	Phospho-FAK (Tyr576) (17X9) Rabbit Monoclonal Antibody
<b>Description</b>	Rabbit Monoclonal Antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB
<b>Reactivity</b>	Human,Mouse,Rat

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Phospho Antibody
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
<b>Purification</b>	Affinity purification

## Immunogen

<b>Gene Name</b>	PTK2
<b>Alternative Names</b>	
<b>Gene ID</b>	5747.0
<b>SwissProt ID</b>	Q05397.A synthetic phosphopeptide corresponding to residues surrounding Tyr576 of human FAK

## Application

<b>Dilution Ratio</b>	WB: 1:2000-1:10000
<b>Molecular Weight</b>	119kDa

**Product Name: Phospho-FAK (Tyr576) (17X9) Rabbit  
Monoclonal Antibody  
Catalog #: AMRe05905**




---

## Background

## Research Area

## Image Data



Western blot detection of FAK (phospho-Tyr576) in A549,U251,C2C12 cell lysates using FAK (phospho-Tyr576) antibody(1:1000 diluted).

## Note

For research use only.