

**Product Name: Phospho-IRF3 (Ser386) (9P16) Rabbit
Monoclonal Antibody
Catalog #: AMRe05928**



Summary

Production Name	Phospho-IRF3 (Ser386) (9P16) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Phospho Antibody
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol.
Purification	Affinity purification

Immunogen

Gene Name	IRF3
Alternative Names	
Gene ID	3661.0
SwissProt ID	Q14653.

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	47 kDa


Background

**Product Name: Phospho-IRF3 (Ser386) (9P16) Rabbit
Monoclonal Antibody
Catalog #: AMRe05928**



Research Area

Image Data



Western blot analysis of Phospho-IRF3 (S386) in MCF-7 lysates using Phospho-IRF3 (Ser386) antibody.

Note

For research use only.