

**Product Name: Phospho-mTOR (Ser2481) (10P6) Rabbit  
Monoclonal Antibody  
Catalog #: AMRe05952**

---

## Summary

<b>Production Name</b>	Phospho-mTOR (Ser2481) (10P6) Rabbit Monoclonal Antibody
<b>Description</b>	Rabbit Monoclonal Antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB
<b>Reactivity</b>	Human,Mouse,Rat

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Phospho Antibody
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol.
<b>Purification</b>	Affinity purification

## Immunogen

<b>Gene Name</b>	MTOR
<b>Alternative Names</b>	
<b>Gene ID</b>	2475.0
<b>SwissProt ID</b>	P42345.

## Application

<b>Dilution Ratio</b>	WB: 1/500-1/1000
<b>Molecular Weight</b>	289 kDa

## Background

---


**Product Name: Phospho-mTOR (Ser2481) (10P6) Rabbit  
Monoclonal Antibody  
Catalog #: AMRe05952**



---

## Research Area

## Image Data



Western blot analysis of Phospho-mTOR (Ser2481) in C6 lysates using Phospho-mTOR (Ser2481) antibody.

## Note

For research use only.