

Product Name: ADFP (18T3) Rabbit Monoclonal Antibody
Catalog #: AMRe06625

Summary

Production Name	ADFP (18T3) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	PLIN2
Alternative Names	ADFP; Adipophilin; Adipose differentiation related protein; ADRP; Perilipin-2; PLIN2;
Gene ID	123.0
SwissProt ID	Q99541.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	48kDa


Product Name: ADFP (18T3) Rabbit Monoclonal Antibody
Catalog #: AMRe06625

Background

May be involved in development and maintenance of adipose tissue. May be involved in development and maintenance of adipose tissue.

Research Area

Image Data



Western blot analysis of ADFP expression in HepG2 cell lysate.

Note

For research use only.