

Product Name: alpha Tubulin (4D19) Rabbit Monoclonal Antibody
Catalog #: AMRe06803

Summary

Production Name	alpha Tubulin (4D19) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	TUBA4A
Alternative Names	
Gene ID	7277.0
SwissProt ID	P68366.A synthetic peptide of human alpha Tubulin

Application

Dilution Ratio	WB: 1:1000-1:10000
Molecular Weight	50kDa


Background

Product Name: alpha Tubulin (4D19) Rabbit Monoclonal Antibody
Catalog #: AMRe06803



Research Area

Image Data



Western blot detection of alpha Tubulin in HeLa,CHO-K1,Rat Brain lysates using alpha Tubulin antibody(1:500 diluted).

Note

For research use only.