

Product Name: Annexin A10 (14O15) Rabbit Monoclonal Antibody
Catalog #: AMRe06915

Summary

Production Name	Annexin A10 (14O15) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	ANXA10
Alternative Names	annexin 14; annexin A10; annexin-14; ANX10; ANX14;
Gene ID	11199.0
SwissProt ID	Q9UJ72.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	37kDa

Product Name: Annexin A10 (14O15) Rabbit Monoclonal Antibody
Catalog #: AMRe06915




Background

Belongs to the annexin family. Contains 4 annexin repeats.

Research Area

Image Data



Western blot analysis of Annexin A10 expression in HeLa cell lysate.

Note

For research use only.