

**Product Name: Apolipoprotein A II (14E13) Rabbit
Monoclonal Antibody
Catalog #: AMRe07050**

Summary

Production Name	Apolipoprotein A II (14E13) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	APOA2
Alternative Names	APO A2; Apo AII; APOA 2; ApoA II; APOA2; APOAII; Apolipoprotein A II;
Gene ID	336.0
SwissProt ID	P02652.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	11kDa


**Product Name: Apolipoprotein A II (14E13) Rabbit
Monoclonal Antibody
Catalog #: AMRe07050**

Background

May stabilize HDL (high density lipoprotein) structure by its association with lipids, and affect the HDL metabolism. May stabilize HDL (high density lipoprotein) structure by its association with lipids, and affect the HDL metabolism.

Research Area

Image Data



Western blot analysis of Apolipoprotein A II expression in Human plasma lysate.

Note

For research use only.