

Product Name: Aquaporin 2 (13E19) Rabbit Monoclonal Antibody
Catalog #: AMRe07079

Summary

| | |
|------------------------|------------------------------------------------|
| Production Name | Aquaporin 2 (13E19) Rabbit Monoclonal Antibody |
| Description | Rabbit Monoclonal Antibody |
| Host | Rabbit |
| Application | WB |
| Reactivity | Human,Mouse,Rat |

Performance

| | |
|---------------------|--------------------------------------------------------------------------------------------------------------|
| Conjugation | Unconjugated |
| Modification | Unmodified |
| Isotype | IgG |
| Clonality | Monoclonal |
| Form | Liquid |
| Storage | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| Buffer | Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA. |
| Purification | Affinity purification |

Immunogen

| | |
|--------------------------|-------------------------------------------------|
| Gene Name | AQP2 |
| Alternative Names | |
| Gene ID | 359.0 |
| SwissProt ID | P41181.A synthetic peptide of human Aquaporin 2 |

Application

| | |
|-------------------------|------------|
| Dilution Ratio | WB: 1:1000 |
| Molecular Weight | 29kDa |


Background

Product Name: Aquaporin 2 (13E19) Rabbit Monoclonal Antibody
Catalog #: AMRe07079



Research Area

Image Data



Western blot detection of Aquaporin 2 in Raoms,Raji cell lysates using Aquaporin 2 antibody(1:1000 diluted).

Note

For research use only.