

Product Name: Beclin 1 (18O10) Rabbit Monoclonal Antibody
Catalog #: AMRe07527



Summary

Production Name	Beclin 1 (18O10) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol.
Purification	Affinity purification

Immunogen

Gene Name	BECN1
Alternative Names	
Gene ID	8678.0
SwissProt ID	Q14457.

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	52 kDa


Background

Product Name: Beclin 1 (18O10) Rabbit Monoclonal Antibody
Catalog #: AMRe07527



Research Area

Image Data



Western blot analysis of Beclin 1 in HEK293, NIH/3T3, C6 lysates using Beclin 1 antibody.

Note

For research use only.