

Product Name: BRCC36 (19W3) Rabbit Monoclonal Antibody
Catalog #: AMRe07645



Summary

Production Name	BRCC36 (19W3) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	BRCC3
Alternative Names	
Gene ID	79184.0
SwissProt ID	P46736.A synthetic peptide of human BRCC36

Application

Dilution Ratio	WB: 1:2000-1:10000
Molecular Weight	36kDa


Background

Product Name: BRCC36 (19W3) Rabbit Monoclonal Antibody
Catalog #: AMRe07645



Research Area

Image Data



Western blot detection of BRCC36 in Rat Brain,C6,3T3,HeLa cell lysates using BRCC36 antibody(1:1000 diluted).

Note

For research use only.