

Product Name: BUBR1 (3N1) Rabbit Monoclonal Antibody
Catalog #: AMRe07697

Summary

Production Name	BUBR1 (3N1) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	BUB1B
Alternative Names	
Gene ID	701.0
SwissProt ID	O60566.A synthetic peptide of human BubR1

Application


Dilution Ratio	WB: 1:1000
Molecular Weight	120kDa

Background

Product Name: BUBR1 (3N1) Rabbit Monoclonal Antibody
Catalog #: AMRe07697

Research Area

Image Data



Western blot detection of BubR1 in HeLa,A549,MCF-7,U251 cell lysates using BubR1 antibody(1:1000 diluted).

Note

For research use only.