

Product Name: C3orf38 (11N12) Rabbit Monoclonal Antibody
Catalog #: AMRe07748



Summary

Production Name	C3orf38 (11N12) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	C3orf38
Alternative Names	C3orf38;
Gene ID	285237.0
SwissProt ID	Q5JPI3.Recombinant protein of human C3orf38

Application

Dilution Ratio	WB: 1:1000-1:5000
Molecular Weight	38kDa

Background


Product Name: C3orf38 (11N12) Rabbit Monoclonal Antibody
Catalog #: AMRe07748



May be involved in apoptosis regulation. May be involved in apoptosis regulation.

Research Area

Image Data



Western blot detection of C3orf38 in K562,HeLa cell lysates using C3orf38 antibody(1:1000 diluted).

Note

For research use only.