

Product Name: Calbindin (12D5) Rabbit Monoclonal Antibody
Catalog #: AMRe07845



Summary

Production Name	Calbindin (12D5) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% New type preservative N and 0.05% BSA
Purification	Affinity purification

Immunogen

Gene Name	Calb1
Alternative Names	
Gene ID	12307.0
SwissProt ID	P12658.

Application

Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100
Molecular Weight	30 kDa


Product Name: Calbindin (12D5) Rabbit Monoclonal Antibody
Catalog #: AMRe07845




Background

Research Area

Image Data



Immunohistochemistry analysis of paraffin-embedded mouse kidney using Calbindin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Calbindin in rat Brain lysates using Calbindin antibody.

Note

For research use only.