

Product Name: Cathelicidin (8H15) Rabbit Monoclonal Antibody
Catalog #: AMRe08007



Summary

Production Name	Cathelicidin (8H15) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	CAMP
Alternative Names	
Gene ID	820.0
SwissProt ID	P49913.A synthetic peptide of human Cathelicidin

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	19kDa


Background

Product Name: Cathelicidin (8H15) Rabbit Monoclonal Antibody
Catalog #: AMRe08007



Research Area

Image Data



Western blot detection of Cathelicidin/CLP in human serum using Cathelicidin/CLP antibody.

Note

For research use only.