

Product Name: CCL23 (3Y14) Rabbit Monoclonal Antibody
Catalog #: AMRe08139

Summary

Production Name	CCL23 (3Y14) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	CCL23
Alternative Names	CCL23; CKB8; Hmrp-2a; MIP3; MPIF1; SCYA23;
Gene ID	6368.0
SwissProt ID	P55773.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	13kDa


Product Name: CCL23 (3Y14) Rabbit Monoclonal Antibody
Catalog #: AMRe08139

Background

Shows chemotactic activity for monocytes, resting T-lymphocytes, and neutrophils, but not for activated lymphocytes. Inhibits proliferation of myeloid progenitor cells in colony formation assays. This protein can bind heparin. Binds CCR1. CCL23(19-99), CCL23(22-99), CCL23(27-99), CCL23(30-99) are more potent chemoattractants than the small-inducible cytokine A23. Shows chemotactic activity for monocytes, resting T- lymphocytes, and neutrophils, but not for activated lymphocytes. Inhibits proliferation of myeloid progenitor cells in colony formation assays. This protein can bind heparin. Binds CCR1. CCL23(19-99), CCL23(22-99), CCL23(27-99), CCL23(30-99) are more potent chemoattractants than the small-inducible cytokine A23.

Research Area

Image Data



Western blot analysis of CCL23 expression in HepG2 cell lysate.

Note

For research use only.