

Product Name: CCL27 (19J3) Rabbit Monoclonal Antibody
Catalog #: AMRe08142

Summary

Production Name	CCL27 (19J3) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	CCL27
Alternative Names	ALP; Ccl27; CTACK; CTAK; ILC; PESKY; SCYA27; Small inducible cytokine A27;
Gene ID	10850.0
SwissProt ID	Q9Y4X3.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	13kDa


Product Name: CCL27 (19J3) Rabbit Monoclonal Antibody
Catalog #: AMRe08142

Background

Chemotactic factor that attracts skin-associated memory T-lymphocytes. May play a role in mediating homing of lymphocytes to cutaneous sites. Binds to CCR10. Chemotactic factor that attracts skin-associated memory T- lymphocytes. May play a role in mediating homing of lymphocytes to cutaneous sites. Binds to CCR10.

Research Area

Image Data



Western blot analysis of CCL27 expression in Human CCL27 recombinant protein cell lysate.

Note

For research use only.