

Product Name: CD226 (7Z13) Rabbit Monoclonal Antibody
Catalog #: AMRe08289

Summary

Production Name	CD226 (7Z13) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	Cd226
Alternative Names	
Gene ID	225825.0
SwissProt ID	Q8K4F0.Recombinant protein of mouse CD226

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	38kDa


Background

Product Name: CD226 (7Z13) Rabbit Monoclonal Antibody
Catalog #: AMRe08289



Research Area

Image Data



Western blot detection of CD226 in mouse spleen cell lysates using CD226 antibody(1:1000 diluted).

Note

For research use only.