

Product Name: CD229 (3L5) Rabbit Monoclonal Antibody
Catalog #: AMRe08290

Summary

Production Name	CD229 (3L5) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	LY9
Alternative Names	
Gene ID	4063.0
SwissProt ID	Q9HBG7.Recombinant protein of human CD229

Application


Dilution Ratio	WB: 1:1000
Molecular Weight	72kDa

Background

Product Name: CD229 (3L5) Rabbit Monoclonal Antibody
Catalog #: AMRe08290

Research Area

Image Data



Western blot detection of CD229 in HeLa,A549,HL-60,U251 cell lysates using CD229 antibody(1:1000 diluted).

Note

For research use only.