

Product Name: CD81 (18O3) Rabbit Monoclonal Antibody
Catalog #: AMRe08464



Summary

Production Name	CD81 (18O3) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	CD81
Alternative Names	
Gene ID	975.0
SwissProt ID	P60033.Recombinant protein of human CD81

Application

Dilution Ratio	WB: 1:1000-1:5000
Molecular Weight	26kDa


Background

Product Name: CD81 (1803) Rabbit Monoclonal Antibody
Catalog #: AMRe08464



Research Area

Image Data



Western blot detection of CD81 in C2C12, 3T3, MCF-7 cell lysates using CD81 antibody(1:1000 diluted).

Note

For research use only.