

Product Name: CD83 (7K14) Rabbit Monoclonal Antibody
Catalog #: AMRe08468

Summary

Production Name	CD83 (7K14) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	CD83
Alternative Names	BL11 PEN; CD83; HB15; hCD83;
Gene ID	9308.0
SwissProt ID	Q01151.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	23kDa


Product Name: CD83 (7K14) Rabbit Monoclonal Antibody
Catalog #: AMRe08468

Background

May play a significant role in antigen presentation or the cellular interactions that follow lymphocyte activation. May play a significant role in antigen presentation or the cellular interactions that follow lymphocyte activation.

Research Area

Image Data



Western blot analysis of CD83 expression in MOLT4 cell lysate.

Note

For research use only.