

Product Name: CD98 (2G13) Rabbit Monoclonal Antibody
Catalog #: AMRe08492

Summary

Production Name	CD98 (2G13) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	SLC3A2
Alternative Names	4T2HC; CD98; CD98HC; MDU1; NACAE; Slc3a2;
Gene ID	6520.0
SwissProt ID	P08195.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	68kDa

Background

Required for the function of light chain amino-acid transporters. Involved in sodium-independent, high-affinity transport of large neutral amino acids such as phenylalanine, tyrosine, leucine, arginine and tryptophan. Involved in guiding and targeting of LAT1 and LAT2 to the plasma membrane. Component of several heterodimeric complexes involved in amino acid transport (PubMed: [11557028](http://www.uniprot.org/citations/11557028)), PubMed: [9829974](http://www.uniprot.org/citations/9829974), PubMed: [9751058](http://www.uniprot.org/citations/9751058), PubMed: [10391915](http://www.uniprot.org/citations/10391915), PubMed: [10574970](http://www.uniprot.org/citations/10574970), PubMed: [11311135](http://www.uniprot.org/citations/11311135), PubMed: [30341327](http://www.uniprot.org/citations/30341327)). The precise substrate specificity depends on the other subunit in the heterodimer (PubMed: [9829974](http://www.uniprot.org/citations/9829974), PubMed: [9751058](http://www.uniprot.org/citations/9751058), PubMed: [10391915](http://www.uniprot.org/citations/10391915), PubMed: [10574970](http://www.uniprot.org/citations/10574970), PubMed: [30867591](http://www.uniprot.org/citations/30867591), PubMed: [10903140](http://www.uniprot.org/citations/10903140)). The complexes function as amino acid exchangers (PubMed: [11557028](http://www.uniprot.org/citations/11557028), PubMed: [10903140](http://www.uniprot.org/citations/10903140), PubMed: [12117417](http://www.uniprot.org/citations/12117417), PubMed: [12225859](http://www.uniprot.org/citations/12225859), PubMed: [30867591](http://www.uniprot.org/citations/30867591)). The homodimer functions as sodium-independent, high-affinity transporter that mediates uptake of large neutral amino acids such as phenylalanine, tyrosine, L-DOPA, leucine, histidine, methionine and tryptophan (PubMed: [9751058](http://www.uniprot.org/citations/9751058), PubMed: [11557028](http://www.uniprot.org/citations/11557028), PubMed: [11311135](http://www.uniprot.org/citations/11311135), PubMed: [11564694](http://www.uniprot.org/citations/11564694), PubMed: [12117417](http://www.uniprot.org/citations/12117417), PubMed: [12225859](http://www.uniprot.org/citations/12225859), PubMed: [25998567](http://www.uniprot.org/citations/25998567), PubMed: [30867591](http://www.uniprot.org/citations/30867591)). The heterodimer formed by SLC3A2 and SLC7A6 or SLC3A2 and SLC7A7 mediates the uptake of dibasic amino acids (PubMed: [9829974](http://www.uniprot.org/citations/9829974), PubMed: [10903140](http://www.uniprot.org/citations/10903140)).


Product Name: CD98 (2G13) Rabbit Monoclonal Antibody
Catalog #: AMRe08492

target="_blank">10903140). The heterodimer with SLC7A5/LAT1 mediates the transport of thyroid hormones triiodothyronine (T3) and thyroxine (T4) across the cell membrane (PubMed:11564694, PubMed:12225859). The heterodimer with SLC7A5/LAT1 is involved in the uptake of toxic methylmercury (MeHg) when administered as the L-cysteine or D,L- homocysteine complexes (PubMed:12117417). The heterodimer with SLC7A5/LAT1 is involved in the uptake of leucine (PubMed:25998567, PubMed:30341327). When associated with LAPTM4B, the heterodimer with SLC7A5/LAT1 is recruited to lysosomes to promote leucine uptake into these organelles, and thereby mediates mTORC1 activation (PubMed:25998567). The heterodimer with SLC7A5/LAT1 may play a role in the transport of L-DOPA across the blood-brain barrier (By similarity). The heterodimer formed by SLC3A2 and SLC7A5/LAT1 or SLC3A2 and SLC7A8/LAT2 is involved in the cellular activity of small molecular weight nitrosothiols, via the stereoselective transport of L-nitrosocysteine (L-CNSO) across the transmembrane (PubMed:15769744). Together with ICAM1, regulates the transport activity of SLC7A8/LAT2 in polarized intestinal cells by generating and delivering intracellular signals (PubMed:12716892). Required for targeting of SLC7A5/LAT1 and SLC7A8/LAT2 to the plasma membrane and for channel activity (PubMed:9751058, PubMed:11311135, PubMed:30867591). Plays a role in nitric oxide synthesis in human umbilical vein endothelial cells (HUVECs) via transport of L-arginine (PubMed:14603368). May mediate blood-to-retina L-leucine transport across the inner blood-retinal barrier (By similarity).

Research Area

Image Data

Product Name: CD98 (2G13) Rabbit Monoclonal Antibody
Catalog #: AMRe08492



Western blot analysis of CD98 expression in HepG2 cell lysate.

Note

For research use only.