

Product Name: Collagen VI (15E2) Rabbit Monoclonal Antibody
Catalog #: AMRe09217



Summary

Production Name	Collagen VI (15E2) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	
Alternative Names	Alpha 1 (VI) chain (61 AA); COL6A1; COL6A2; COL6A3; Collagen alpha 2(VI) chain; Collagen alpha 3(VI) chain; OPLL; PP3610;
Gene ID	
SwissProt ID	P12109/P12110/P12111.

Application

Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	147kDa

Product Name: Collagen VI (15E2) Rabbit Monoclonal Antibody
Catalog #: AMRe09217




Background

Collagen VI acts as a cell-binding protein.

Research Area

Image Data



Western blot analysis of Collagen VI expression in Human brain lysate.

Note

For research use only.