

Product Name: CRYAB (9C10) Rabbit Monoclonal Antibody
Catalog #: AMRe09436

Summary

Production Name	CRYAB (9C10) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	CRYAB
Alternative Names	AACRYA; Alpha B crystallin; Alpha-crystallin B chain; CRYA2; Cryab; Crystallin alpha B; CTPP2; Heat shock 20 kD like protein; HspB5;
Gene ID	1410.0
SwissProt ID	P02511.A synthetic peptide of human Alpha B Crystallin

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	20kDa


Product Name: CRYAB (9C10) Rabbit Monoclonal Antibody
Catalog #: AMRe09436

Background

May contribute to the transparency and refractive index of the lens. May contribute to the transparency and refractive index of the lens. Has chaperone-like activity, preventing aggregation of various proteins under a wide range of stress conditions.

Research Area

Image Data



Western blot detection of Alpha B Crystallin in C6,HeLa cell lysates using Alpha B Crystallin antibody(1:1000 diluted).

Note

For research use only.